

Introduction

To promote re-epithelialisation, a wound must remain moist however, if a dressing retains too much moisture, maceration of the surrounding healthy tissue can occur (Junker *et al.* 2013). Wound dressings should have sufficient permeability to moisture vapour to prevent backward migration collecting within the dressing (Zehrer *et al.* 2014) however, entrapment of liquid can lead to serious health consequences for skin integrity. The Moisture Vapour Transmission Rate (MVTR) test method measures both water vapour and water contact for wound dressings using mass differential.

Methodology

Three wound dressings were selected for MVTR assessment. The dressings were preconditioned at $21^{\circ}\text{C} \pm 2^{\circ}\text{C}$ at a relative humidity (RH) of $60 \pm 15\% \text{RH}$ for 48 hours. Dressings were clamped to the cylinders and filled with distilled water, leaving a gap of 5 ± 1 mm between the water and the dressing. To assess the dressing in contact with water vapour, each device was weighed and incubated at $37 \pm 1^{\circ}\text{C}$ for 18 - 24 hours. To test the MVTR in contact with water, the Paddington cups were inverted and incubated at $37 \pm 1^{\circ}\text{C}$ for 4 hours.




Figure 1. The average moisture vapour transmission rate of three test dressings in contact with water (4 hours contact time) and in contact with vapour (24 hours contact time).



Figure 2. Cylinders used for the MVTR test method

Results

The observed ranges for dressings A, B and C when in contact with water were $40185.84 \pm 12926.92 \text{ gm}^{-2}24\text{h}^{-1}$, $2342.52 \pm 2409.61 \text{ gm}^{-2}24\text{h}^{-1}$ and $9740.40 \pm 6886.89 \text{ gm}^{-2}24\text{h}^{-1}$, respectively. The control dressing in contact with water resulted in a range of $27.24 \pm 47.36 \text{ gm}^{-2}24\text{h}^{-1}$.

The observed ranges for dressings A, B and C when in contact with water vapour were $926.40 \pm 1956.67 \text{ gm}^{-2}24\text{h}^{-1}$, $564.96 \pm 202.73 \text{ gm}^{-2}24\text{h}^{-1}$ and $563.50 \pm 380.19 \text{ gm}^{-2}24\text{h}^{-1}$, respectively. The control dressing when in contact with water vapour resulted in a range of $469.18 \pm 315.52 \text{ gm}^{-2}24\text{h}^{-1}$.

Discussion and Conclusions

According to Xu *et al.* (2016), an MVTR of $2028.30 \pm 237.8 \text{ gm}^{-2}24\text{h}^{-1}$ equates to the optimal moist environment which promotes wound healing. Dressing A demonstrated the highest range for both water vapour and water contact which was expected as this test sample is the most porous of the samples tested and would not be used for an exuding wound in the real world. Dressing B demonstrated the lowest range for both water vapour and water contact which was expected as this is test sample is categorized as a permeable film dressing and would be used for an exuding wound in the real world. Dressing C demonstrated a higher MVTR result compared to dressing B for both water contact and water vapour which was expected due to being suited to helping manage moderate to high exuding wounds. According to Xu *et al.* (2016), an MVTR of $2028.30 \pm 237.8 \text{ gm}^{-2}24\text{h}^{-1}$ equates to the optimal moist environment which promotes wound healing.

References

- Junker, J.P., Kamel, R.A., Caterson, E.J. and Eriksson, E., 2013. Clinical impact upon wound healing and inflammation in moist, wet, and dry environments. *Advances in wound care*, 2(7), pp.348-356.
- Zehrer, C.L., Holm, D., Solfest, S.E. and Walters, S.A., 2014. A comparison of the in vitro moisture vapour transmission rate and in vivo fluid-handling capacity of six adhesive foam dressings to a newly reformulated adhesive foam dressing. *International wound journal*, 11(6), pp.681-690.
- Xu, R., Xia, H., He, W., Li, Z., Zhao, J., Liu, B., Wang, Y., Lei, Q., Kong, Y., Bai, Y. and Yao, Z., 2016. Controlled water vapor transmission rate promotes wound-healing via wound re-epithelialization and contraction enhancement. *Scientific reports*, 6, p.24596.

This project was carried out by Perfectus Biomed Group

Perfectus Biomed Group, Sci-Tech Daresbury, Keckwick Lane, Cheshire, UK, WA4 4AD Tel +44 1925 864838

Perfectus Biomed Group, 3545 South Park Drive, Jackson, Wyoming 83001, USA Tel: +1 307 264 1572

Email: info@perfectusbiomed.com