

Introduction

Over 65% of bacterial infections treated clinically in the developed world are known to be caused by organisms growing in biofilms. Biofilms are matrix-enclosed communities of bacteria which are significantly more resistant to host defences and to conventional therapies compared to their planktonic counterparts. Periprosthetic infections can lack a definitive diagnosis because matrix-enclosed sessile bacteria are less immunogenic and elicit a reduced inflammatory response compared to planktonic cells. Implantable devices are highly susceptible to infection and biofilm formation. Once a biofilm has been established in a periprosthetic joint can be difficult to diagnose and eradicate. Successful treatment of periprosthetic joint infection requires optimal surgical procedures combined with long-term antimicrobial therapy directed against surface-adhering microorganisms. Given the limited efficacy of traditional antibiotics in implant-associated infections, novel strategies such as releasing antibiotics directly to the site of infection are a promising future option for biofilm prevention and eradication.

Aim

To assess the ability of Stimulan® beads containing a mixture of Vancomycin and Gentamicin or Vancomycin and Tobramycin to release a combination of antibiotics and effectively eradicate pre-formed biofilms *in vitro*.

Method

- Single species *Pseudomonas aeruginosa* and *Staphylococcus aureus* biofilms were established on polycarbonate coupons within a CDC biofilm reactor for 72 hours (Figure 1).
- Established biofilms were removed from the reactor and rinsed in phosphate buffered saline to remove planktonic organisms.
- Washed biofilms were exposed to a challenge plate containing suspended Stimulan® beads enclosing a mixture of Vancomycin and Gentamicin (VG) or Vancomycin and Tobramycin (VT) for 24 hours at 37°C ± 2°C.
- Positive and negative controls were tested concurrently.
- Following exposure, remaining attached organisms were recovered by sonication. Serial dilutions and plate counts were performed on the resultant suspensions. Log reductions compared to the negative control were reported.




Figure 1. Photographs of (A) CDC reactor, (B) Coupon holding rods, (C) CDC reactor containing 72-hour pre-formed single species biofilms.

Results

An average of $6.60 \pm 0.23 \text{ Log}_{10} \text{CFU mL}^{-1}$ viable *S. aureus* and $6.78 \pm 0.12 \text{ Log}_{10} \text{CFU mL}^{-1}$ viable *P. aeruginosa* were recovered from the negative controls (Figure 2). No viable *S. aureus* were recovered from biofilms treated with the positive control or Stimulan® beads containing a mixture of Vancomycin and Gentamicin or Vancomycin and Tobramycin. No viable *P. aeruginosa* were recovered from biofilms treated with the positive control or Stimulan® beads containing a mixture of Vancomycin and Gentamicin. An average of $0.33 \pm 0.52 \text{ Log}_{10} \text{CFU mL}^{-1}$ (below 1 Log) viable *P. aeruginosa* were recovered from biofilms following treatment with Stimulan® beads containing a mixture of Vancomycin and Tobramycin.




Figure 2. Quantity of total viable microorganisms recovered following 24-hour treatment of 72-hour pre-formed *Staphylococcus aureus* (A) and *Pseudomonas aeruginosa* (B) biofilms.

Discussion

In vitro biofilm test methods are more representative of the real-world scenario compared to their planktonic equivalents. Biofilm microorganisms are attached to a surface, have developed under a dynamic environment and reside within extra polymeric substances. Exposure of the biofilm to Stimulan® beads containing a mixture of Vancomycin and Gentamicin or Vancomycin and Tobramycin resulted no viable organisms being recovered from pre-formed biofilms in the test method described. Further testing is required to confirm clinical performance. Eradication of biofilm infections may result in improved patient outcomes and a decrease in healthcare costs. Demonstration of biofilm eradication *in vitro* provides a strong basis for future work to show the potential for effective treatment clinically.

References:

1. Training the Biofilm Generation – a Tribute to J.W. Costerton. McClean RJC, Lam, JS and Graham LL. *J. Bacteriol.* 2012 Dec; 194(24): 6706–6711.
2. Bacterial biofilms: a common cause of persistent infections. Costerton JW, Stewart PS, Greenberg EP. *Science.* 1999 May 21; 284(5418):1318-22.
3. Biofilm theory can guide the treatment of device-related orthopaedic infections. Costerton JW. *Clin Orthop Relat Res.* 2005 Aug; (437):7-11.